

7" Love-Gruenwald IVD w/screw, angled up

SKU: 70-1124

EXAMPLE



Style	IVD
Secondary Style	Love-Gruenwald
Additional Style	screw, teeth
Bite Width	3mm
Bite Length	10mm
Angle	up
Working Length	7"
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

Love-Gruenwald IVD, screw, teeth, up, 3x10mm, 7"

Picture of 70-1054 displayed as an example.