

8" Tarpon Kerrison Rongeur 2mm Pin Ejector 15mm bite Stainless Steel

SKU: 70-1862



Style	Kerrison
Secondary Style	Tarpon
Bite Width	2mm
Bite Length	15mm
Angle	40 up
Working Length	8"
Handle Style	Midsized
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

8" Tarpon Kerrison Rongeur 2mm Pin Ejector 15mm bite Stainless Steel is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.