

8" Micro Kerrison Rongeur, 1mm Super-Coat

SKU: 70-0162SC



Style	Kerrison
Secondary Style	Micro
Bite Width	1mm
Bite Length	9mm
Angle	40 up
Working Length	8"
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

8" Micro Kerrison Rongeur, 1mm Super-Coat is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.