

6" Bayonet Kerrison Rongeur 3mm Pin Ejector Micro Handle Panther

SKU: 70-2281P



Style	Kerrison
Secondary Style	Micro
Additional Style	Bayonet
Bite Width	3mm
Bite Length	11mm
Angle	40 up
Working Length	6"
Handle Style	Small
Material	Stainless Steel
Single Use or Reusable	Reusable

6" Bayonet Kerrison Rongeur 3mm Pin Ejector Micro Handle Panther is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.