

6" Bayonet Kerrison Rongeur 1mm Spurling Handle Panther

SKU: 70-2320P



Style	Kerrison
Secondary Style	Spurling
Additional Style	Bayonet
Bite Width	1mm
Bite Length	15mm
Angle	40 up
Working Length	6"
Handle Style	Large
Material	Stainless Steel
Single Use or Reusable	Reusable

6" Bayonet Kerrison Rongeur 1mm Spurling Handle Panther is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.