

Single-Use Long Frazier Suction 15fr 20cm Polished Tip

SKU: S265



Diameter	15fr
Working Length	20cm
Style	Frazier
Tip Style	Polished
Control Hole Style	No Control Hole
Bend Angle (deg)	30
Addl Style	Long
Material	Stainless Steel Tube/Plastic Handle
Single Use or Reusable	Single-Use
Unit of Sale	BX(25)

Long Suction 15fr 20cm No Control Hole Single Use / 25 per box

Picture of S021 displayed as an example.