

Spinal Fusion Curette - size 2-0

SKU: 70-7251

EXAMPLE



Style	Brun
Secondary Style	Spinal Fusion
Tip Size	2-0
Tip Type	Oval
Overall Length	9"
Sharp Blunt	Blunt
Curve Angle	Straight
Malleable Rigid	Rigid
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

Spinal Fusion Curette, hex handle,str,9",sz 2-0

Picture of 70-7253 displayed as an example.