

Kerrison Rongeur, Micro, 9" (225mm) 3mm, 90° up

Stainless Steel

SKU: 70-0344

Style	Kerrison
Secondary Style	Micro
Bite Width	3mm
Bite Length	9mm
Angle	90 up
Working Length	9"
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

Kerrison Rongeur, Micro, 9" (225mm) 3mm, 90° up Stainless Steel

is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.