

Kerrison Rongeur, Micro, 9" (225mm) 5mm, 40° up Panther

SKU: 70-0358P

Style	Kerrison
Secondary Style	Micro
Bite Width	5mm
Bite Length	9mm
Angle	40 up
Working Length	9"
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

Kerrison Rongeur, Micro, 9" (225mm) 5mm, 40° up Panther is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.