

Kerrison Rongeur, Tarpon, 9" (225mm) 3mm, 40° up Panther

SKU: 70-1953P



Style	Kerrison
Secondary Style	Tarpon
Bite Width	3mm
Bite Length	15mm
Angle	40 up
Working Length	9"
Handle Style	Midsized
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

Kerrison Rongeur, Tarpon, 9" (225mm) 3mm, 40° up Panther is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.