

Kerrison Rongeur, Tarpon, 10" (250mm) 5mm, 40° up, Pin Ejector Stainless Steel

SKU: 70-1976

Style	Kerrison
Secondary Style	Tarpon
Bite Width	5mm
Bite Length	15mm
Angle	40 up
Working Length	10"
Handle Style	Midsize
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

Kerrison Rongeur, Tarpon, 10" (250mm) 5mm, 40° up, Pin Ejector Stainless Steel is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.