

Kerrison Rongeur, Bayonet, 6" (150mm) 6mm, 90° up, Pin Ejector Spurling Handle, Panther

SKU: 70-2358P

Style	Kerrison
Secondary Style	Spurling
Additional Style	Bayonet
Bite Width	6mm
Bite Length	15mm
Angle	90 up
Working Length	6"
Handle Style	Large
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)

Kerrison Rongeur, Bayonet, 6" (150mm) 6mm, 90° up, Pin Ejector Spurling Handle, Panther is part of our premium collection of German-crafted Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.