

Rhoton Set case and 19 instruments titanium

SKU: 71-5123T



Material	Titanium
Single Use or Reusable	Reusable

Set includes:

SKU	Quantity	Description
71-5101T	1	Round Diss, 1mm, 19cm, ang, titanium
71-5102T	1	Round Diss, 2mm, 19cm, titanium
71-5103T	1	Round Diss, 3mm, 19cm, titanium
71-5104T	1	Elevator, cvd, titanium
71-5105T	1	Elevator, ang, titanium Rhoton
71-5106T	1	Spatula Diss, Titanium sm

71-5107T	1	Spatula Diss, titanium, med
71-5108T	1	Spatula Diss, titanium, lg
71-5109T	1	Rt. Angle Hook, semi shrp, titanium
71-5110T	1	Rt. Angle Hook, Titanium blnt, titanium
71-5111T	1	Needle Diss, semi-sharp, 45deg ang, titanium
71-5112T	1	Needle Diss, semi-sharp, str, titanium
71-5113T	1	Micro Curette, str, 1x2mm cup, titanium
71-5114T	1	Micro Curette, ang, 1x2mm cup, titanium
71-5115T	1	Ball Diss, 3mm, 90deg ang, titanium
71-5116T	1	Rhoton Ball Dissector, 5mm, 90 deg, titanium
71-5117T	1	Rhoton Ball Dissector, 5mm, 45 deg, titanium
71-5118T	1	Rhoton Ball Dissector, 8mm, 45 deg, titanium
71-5119T	1	Rhoton Ball Dissector, str, titanium
71-5122	1	Rhoton Case for 19 instruments 10" x 10 3/4" x 1" deep