

CleanWave Kerrison, 6mm, 40 deg up, Panther

SKU: 70-1866CWP

EXAMPLE



Style	Kerrison
Additional Style	Thin Footplate
Ejector	Ejector
Bite Width	6mm
Bite Length	15mm
Angle	40 up
Length	8" (200mm)
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)
Performance Coating	Panther

CleanWave Kerrison, 6mm, 40 deg up, Panther is part of our premium collection of Kerrison rongeurs, designed for consistent performance. We stand behind our instruments with a quality guarantee on the workmanship and materials for the life of each instrument.

Picture of 70-1865CWP displayed as an example.