104 Sommerfield Drive, Gordonsville, VA 22942
800.210.2677 | bossinstruments.com

Sinus Curette, 2x5mm cup, curved 90°, 7"

SKU: 93-5548



Style	Sinus
Tip Size	2mm x 5mm
Tip Type	Oval
Overall Length	7.5"
Sharp | Blunt	Blunt
Curve | Angle	Curved 90
Malleable | Rigid	Rigid
Material	Stainless Steel
Single Use or Reusable	Reusable
Unit of Sale	EA(1)